




# LEO

## HANDHELD ANALYZER FOR IVF APPLICATIONS




## QUICK INSTRUCTIONS

### USER CALIBRATION WITH GAS CYLINDER




### GAS CYLINDER WITH PRESSURE REGULATOR


- 1 Connect the **clear connector (a)** of Tube-A to LEO gas inlet port labeled with icon
  - 2 Put the cap on LEO gas inlet port labeled with icon
  - 3 Connect the **red connector (b)** of Tube-A to the **clear connector (c)** of the Calibration Kit
  - 4 Connect the end of Calibration Kit to the Gas Cylinder
- Use the fitting tubes and/or additional tube in order to allow the compatibility of the supplied tubing with the gas port of your gas cylinders.
- 5 Use the supplied female luer to connect a purge tube to LEO gas outlet port. Place the purge tube in a well ventilated environment
  - 6 Set 1.0 barg on the pressure regulator of the Gas Cylinder
- Do not exceed 1.5 barg
- 7 Start the calibration (refer to LEO User Manual, paragraph 13)




### GAS CYLINDER WITH FLOW REGULATOR




- 1 Connect the **clear connector (a)** of Tube-A to LEO gas inlet port labeled with icon
  - 2 Put the cap on LEO gas inlet port labeled with icon
  - 3 Connect the **red connector (b)** of Tube-A to the **green connector (c)** of Tube-B
  - 4 Connect the **end of Tube-B** to the Calibration Gas Regulator
- Use the fitting tubes and/or additional tube in order to allow the compatibility of the supplied tubing with the gas port of your gas cylinders.
- 5 Use the supplied female luer to connect a purge tube to LEO gas outlet port. Place the purge tube in a well ventilated environment
  - 6 Open the Calibration Gas Regulator by using the valve
- Set the flow rate value at maximum 0.4 l/min.
- 7 Start the calibration (refer to LEO User Manual, paragraph 13)




### USER CALIBRATION WITH GAS AT ATMOSPHERIC PRESSURE



- 1 Connect the **clear connector (a)** of Tube-A to LEO gas inlet port labeled with icon
  - 2 Put the cap on LEO gas inlet port labeled with icon
  - 3 Connect the **red connector (b)** of Tube-A to the **green connector (c)** of Tube-B
  - 4 Connect the **end of Tube-B** to the incubator or the point where measurement is required
- Use the fitting tubes and/or additional tube in order to allow the compatibility of the supplied tubing with the gas port of your incubators.
- 5 Use the supplied female luer to connect a purge tube to LEO gas outlet port. Place the purge tube in a well ventilated environment
- 6 Start the calibration (refer to LEO User Manual, paragraph 13)



### GAS MEASUREMENT IN DIFFUSION MODE




- 1 Connect the **clear connector (a)** of Tube-A to LEO gas inlet port labeled with icon
- 2 Put the cap on LEO gas inlet port labeled with icon
- 3 Connect the **red connector (b)** of Tube-A to the **green connector (c)** of Tube-B
- 4 Connect the **end of Tube-B** to the incubator or the point where measurement is required

Use the fitting tubes and/or additional tube in order to allow the compatibility of the supplied tubing with the gas port of your incubators


Do not exceed 300 mbarg upstream LEO

- 5 Use the supplied female luer to connect a purge tube to LEO gas outlet port. Place the purge tube in a well ventilated environment
- 6 Start a Gas Single point (refer to LEO User manual paragraph 10.1.3) or a Gas Logging (refer to LEO User Manual, paragraph 10.1.4)



Do not connect purge tube without luer connector


### GAS MEASUREMENT IN ASPIRATION MODE



- 1 Connect the **clear connector (a)** of Tube-A to LEO gas inlet port labeled with icon
- 2 Put the cap on LEO gas inlet port labeled with icon
- 3 Connect the **red connector (b)** of Tube-A to the **green connector (c)** of Tube-B
- 4 Connect the **end of Tube-B** to the incubator or the point where measurement is required


Use the fitting tubes and/or additional tube in order to allow the compatibility of the supplied tubing with the gas port of your incubators.

- 5 Use the supplied female luer to connect a purge tube to LEO gas outlet port. Place the purge tube in a well ventilated environment
- 6 Start a Gas Single point (refer to LEO User manual paragraph 10.1.3) or a Gas Logging (refer to LEO User Manual, paragraph 10.1.4)



Do not connect purge tube without luer connector

### GAS MEASUREMENT WITH MOISTURE TRAP



If the sampling gas is humid, connect the **red connector (b)** of Tube-A and the **green connector (c)** of Tube-B to Moisture Trap

Before starting a Gas Logging, unscrew the condensation trap to make sure the condensation trap top has been removed




The condensation trap top allows reducing the sampling gas time during a Gas Single Point measurement



Replace the condensation trap in its housing.

### CO2 ZERO RESET WITH N2




Do not use the Calibration Kit without Tube-A


Do not connect purge tube without luer connector

### GAS CYLINDER WITH PRESSURE REGULATOR

- 1 Connect the **clear connector (a)** of Tube-A to LEO gas inlet port labeled with icon
  - 2 Put the cap on LEO gas inlet port labeled with icon
  - 3 Connect the **clear connector (c)** of the Calibration Kit to the **red connector (b)** of Tube-A
  - 4 Connect the **end of Calibration Kit** to the Gas Cylinder
- Use the fitting tubes and/or additional tube in order to allow the compatibility of the supplied tubing with the gas port of your gas cylinders.
- 5 Use the supplied female luer to connect a purge tube to LEO gas outlet port. Place the purge tube in a well ventilated environment
  - 6 Set 1.0 barg on the pressure regulator of the pure N2 tank
- Do not exceed 1.5 barg
- 7 Start the CO2 Zero Reset (refer to LEO User Manual, paragraph 13.4.2)




### GAS CYLINDER WITH FLOW REGULATOR




Do not connect purge tube without luer connector


- 1 Connect the **clear connector (a)** of Tube-A to LEO gas inlet port labeled with icon
  - 2 Put the cap on LEO gas inlet port labeled with icon
  - 3 Connect the **red connector (b)** of Tube-A to the **green connector (c)** of Tube-B
  - 4 Connect the **end of Tube-B** to the Calibration Gas Regulator
- Use the fitting tubes and/or additional tube in order to allow the compatibility of the supplied tubing with the gas port of your gas cylinders.
- 5 Use the supplied female luer to connect a purge tube to LEO gas outlet port. Place the purge tube in a well ventilated environment
  - 6 Open the Calibration Gas Regulator by using the valve
- Set the flow rate value at maximum 0.4 l/min.
- 7 Start the CO2 Zero Reset (refer to LEO User Manual, paragraph 13.4.2)



### CO2 ZERO RESET WITH AIR



- 1 Connect the **clear connector (a)** of Tube-A to LEO gas inlet port labeled with icon
- 2 Put the cap on LEO gas inlet port labeled with icon
- 3 Connect the **red connector (b)** of Tube-A to the **clear connector (c)** of Soda Lime Kit
- 4 Use the supplied female luer to connect a purge tube to LEO gas outlet port. Place the purge tube in a well ventilated environment
- 5 Start the CO2 Zero Reset (refer to LEO User Manual, paragraph 13.4.1)



Do not use the Soda Lime Kit without Tube-A

Do not connect purge tube without luer connector